

Compliance Chart

Inflation Pressure [bar]	Expanded Stent Outer Diameter [mm]						
	Ø 2.50	Ø 2.75	Ø 3.00	Ø 3.50	Ø 4.00	Ø 4.50	Ø 5.00
7						4.03	4.48
8	2.26	2.49	2.68	3.05	3.52	4.20	4.68
9	2.34	2.58	2.80	3.20	3.69	4.37	4.85
10	2.42	2.67	2.91	3.35	3.86	4.50	5.00
11	2.50	2.75	3.00	3.50	4.00	4.61	5.12
12	2.57	2.82	3.07	3.62	4.13	4.70	5.23
13	2.63	2.88	3.14	3.71	4.23	4.77	5.31
14	2.68	2.93	3.20	3.80	4.31	4.83	5.37
15	2.73	2.98	3.26	3.86	4.39		
16	2.78	3.02	3.32	3.92	4.46		

NP

Nominal Pressure

RBP

Rated Burst Pressure